

Packed Column GC

Custom Packed GC Columns



Packed GC Column

- Custom-packed to your specifications
 - Quality supports and stationary phases
 - Includes standard brass fittings
 - Fast delivery - made in our own facility
- For Glass GC Columns, see pages 30 - 33*



We make custom-packed GC columns using high-quality materials. Any stationary phase can be chosen from our listings on page 36 - 37, in combination with nearly all listed supports. The use of other stationary phases and those not listed in the Packing Surcharge Table are subject to special quotation.

Columns are supplied as 6" diameter coils, unless instrument make and model is specified or dimensions are specified.


Stationary Phases

Usually the solid support will be coated with a liquid stationary phase. Exceptions are molecular sieves and porous polymers such as HayeSep[®], Porapak[®] or Chromosorb Century Series[®], which are normally used un-coated. For a complete listing of available stationary phases, see pages 36 - 37.

Length of injector arm.



Length of detector arm.



Span, the distance between arms, tube center to tube center.

Solid Supports

Diatomaceous Earth Supports

- Chromosorb[®] G (AW, NAW, AW-DMCS)
- Chromosorb P (AW, NAW, AW-DMCS, HMDS)
- Chromosorb W (AW, NAW, AW-DMCS, HMDS, HP)

Porous Polymer Supports

- Chromosorb 101, 102, 103, 104, 105, 106, 107, 108
- HaySep A, B, C, D, DB, DIP, N, P, Q, R, S, T
- Tenax[®]
- Porapak N, PS, Q, QS, R, S, T

Molecular Sieves

- 4A, 5A, 13X

Miscellaneous Supports

- Chromosorb T (PTFE)
- Carbon Molecular Sieves

Note*

- NAW - not acid washed
- AW - acid washed
- AW-DMCS - acid washed and silanized
- HMDS - silanized with hexamethyl disilazane
- High Performance (HP) - specially treated for superior inertness and column efficiency

**For complete descriptions, see individual listings on pages 34 - 39. For other supports, please call 800-327-3800 for a quote.*

Custom Packed GC Columns

To Order a Custom Column

Use the catalog number from the Item Number Guide to indicate column diameter and approximate length. Then specify the exact length needed, as well as other requirements for the solid support and stationary phase, using the listings below.

Item Number Guide

Length Meters	Length Feet	1/8" OD	1/4" OD	1/16" OD	3/16" OD
0.3 - 1.1	1 - 3.9	455501	455521	455541	455502
1.2 - 2.1	4 - 6.9	455503	455523	455543	455504
2.2 - 3.0	7 - 9.9	455505	455525	455545	455506
3.1 - 3.9	10 - 12.9	455507	455527	455547	455508
4.0 - 4.8	13 - 15.9	455509	455529	455549	455510
4.9 - 5.7	16 - 18.9	455511	455531	455542	455512
5.8 - 6.6	19 - 21.9	455513	455533	455544	455514
6.7 - 7.5	22 - 24.9	455515	455535	455546	455516
7.6 - 8.5	25 - 28	455517	455537	455548	455518

Custom Column Ordering

To completely specify a custom column, please provide the information requested below.

- Item Number
- Length - feet or meters
- Tubing Material - stainless steel
- Stationary Phase** - see listings page 36 - 37. *Phases subject to surcharges; see Table 1*
 - Stationary Phase 1, plus % loading
 - Stationary Phase 2* (if dual phase), plus % loading
- Solid Support - specify material and grade (for example, Chromosorb® W-AW; HayeSep® D)

Expensive supports have surcharges; see Table 2

 - Support
 - Mesh Size - specify one 60/80, 80/100 or 100/120
- Instrument Manufacturer
- Model Number
- Fittings - brass or stainless steel*
- Configuration (optional) - for instruments or shapes not shown on pages 31 - 33, please include a drawing with accurate dimensions

Table 1

Stationary Phase Surcharges

Surcharges apply to the following expensive phases:

- EGSS-X
- OV-225
- OV-275
- OV-330
- OV-1701
- Polyphenylether
- Silars
- Dextsils

Table 2

Support Surcharges

- Chromosorb T
- Carbon Molecular Sieves
- HayeSep
- Porapak®
- Chromosorb 101 - 109
- Tenax®
- Carbopacks

*Surcharge will be added

**Not required with solid adsorbents such as HayeSep, Molecular Sieves, etc.

Packed Column GC

Glass Columns - Silane Treated



Glass Column

- Custom Packed to your specifications
 - Silanized for inertness
 - High quality
 - Available for most instruments
- *For Stainless Steel columns see pages 28 - 29*

CRS's glass columns are made from 1/4" OD precision-bore borosilicate tubes (5 mm OD for Shimadzu instruments), which are silanized for inertness. They are individually packaged for safe shipment. All glass columns include 1/4" brass fittings and graphite ferrules (except columns for Shimadzu instruments, provided with graphite ferrules only).

Column Configuration (pages 31 - 33). The length of the "arms" is indicated by dimensions X and Y. These dimensions along with the span between the arms allows confirmation that the column will fit correctly in your instrument.

Column Installation All of our glass columns are supplied with nuts and graphite ferrules. We recommend graphite ferrules because they are soft, to seal securely around the column without excessive tightening. They will not stick to the glass, for easy removal. Complete installation instructions are included with each column.

Recommended Accessories:

- 1/4" Graphite ferrules (10/pk) - Item 211400

To Order a Packed Column

Specify the "empty" catalog number with a "P" added to the end, and specify the packing (example: 3% OV-225 on Chromosorb® W-HP, 80/100 mesh). Some expensive packing materials are subject to a surcharge; see the Surcharge Table.

1. Catalog Number
2. Stationary Phase (see listings on page 36 - 37)
Expensive phases have surcharges; see Surcharge Table
 - A) Stationary Phase 1
% Loading
 - B) Stationary Phase 2 - Surcharge for dual phase columns
% Loading
3. Solid Support - Specify material and grade (example: Chromosorb W-AW; HayeSep® D)
Expensive supports have surcharges; see Surcharge Table
 - C) Support
 - D) Mesh Size- specify one: 60/80, 80/100, or 100/120

Surcharge Table: Surcharges apply to the following expensive phases:

Stationary Phase Surcharges

- Dexsils
- EGSS-X
- OV-225
- OV-275
- OV-330
- OV-1701
- Polyphenylether
- Silars

Support Surcharges




- Chromosorb T
- Carbon Molecular Sieves
- HayeSep
- Porapak®
- Chromosorb 101 - 109
- Tenax®
- Carbopacks

Glass Columns- Silane Treated

When Ordering a Packed Column

Add a "P" to the empty column part number, for example 220-022P, and specify the packing (example: 3% OV-225 on Chromosorb W-HP, 80/100 mesh). Please provide information on the packing support and stationary phase needed, see page 30.

For Agilent Instruments

Agilent Instrument	Drawing (not to scale)	Column Dimensions	Length	2mm ID Item No.	4mm ID Item No.
5880, 5890, 5987, 6890 Configuration A OCI		X= 280 mm	2ft	220-022	
		Y= 230 mm	3ft	220-032	
		S= 229 mm	4ft	220-042	220-044
			6ft	220-062	220-064
			2m	220-2.0M2	220-2.0M4
			8ft	220-082	
			10ft	220-102	
Configuration B OCI		X= 280 mm	2ft	221-022	221-024
		Y= 180 mm	6ft	221-062	221-064
		S= 229 mm	other	221-S2	221-S4
5830, 5840, 5985 OCI 152 mm span		X= 272 mm	3ft	216-032-6	216-034-6
		Y= 180 mm	4ft	216-042-6	216-044-6
		S= 152 mm	6ft	216-062-6	216-064-6
			8ft	216-082-6	
			10ft	216-102-6	
			other	216-S2-6	216-S4-6
OCI, 229 mm span		X= 272 mm	3ft	216-032-9	
		Y= 182 mm	6ft	216-062-9	216-064-9
		S= 229 mm	8ft	216-082-9	
			10ft	216-102-9	
			other	216-S2-9	216-S4-9
5830, 5840, 5985 NOCI		X= 227 mm	2ft	217-022-6	
		Y= 235 mm	3ft	217-032-6	
		S= 152 mm	6ft	217-062-6	217-064-6
			10ft	217-102-6	
			other	217-S2-6	217-S4-6
NOCI, 229 mm span		X= 227 mm	2ft	217-022-9	
		Y= 235 mm	6ft	217-062-9	
		S= 229 mm	other	217-S2-9	217-S4-9



X = Injector arm
Y = Detector arm
S = Span (distance between arms)

Note
NOCI = not on-column injection
HOCI = heated on-column injection
OCI = on-column injection


Packed Column GC

Glass Columns - Silane Treated

Length of injector arm.



Length of detector arm.



Span, the distance between arms, tube center to tube center.

When Ordering a Packed Column

Add a "P" to the empty column part number, for example 220-022P, and specify the packing (example: 3% OV-225 on Chromosorb W-HP, 80/100 mesh). Please provide information on the packing support and stationary phase needed, see page 30.

For PerkinElmer Instruments



PerkinElmer Instrument	Drawing (not to scale)	Column Dimensions	Length	2mm ID Item No.	4mm ID Item No.
115, 900, 990, F30, 2000 2100, 3920, Sigma Series NOCI		X= 219 mm	2ft	222-022	222-024
		Y= 219 mm	3ft	222-032	222-034
		S= 222 mm	4ft	222-042	222-044
			6ft	222-062	222-064
			8ft	222-082	
			other	222-S2	222-S4
3920, 910, 990, HOCl All glass system		X= 260 mm	3ft	224-032	
		Y= 117 mm	6ft	224-062	
		S= 222 mm	10ft	224-102	
			other	224-S2	224-S4
Sigma Series HOCl		X= 320 mm	3ft	223-032	
		Y= 181 mm	6ft	223-062	223-064
		S= 222 mm	10ft	223-102	
			other	223-S2	223-S4
Sigma 2000		X= 270 mm	3ft	290-032	
		Y= 139 mm	6ft	290-062	290-064
		S= 222 mm	other	290-S2	290-S4
Autosystem 8000, 8300 Series NOCI		X= 173 mm	3ft	226-032	226-034
		Y= 173 mm	6ft	226-062	226-064
		S= 165 mm	other	226-S2	226-S4
Autosystem 8000, 8300 Series OCI		X= 305 mm	3ft	227-032	227-034
		Y= 173 mm	6ft	227-062	227-064
		S= 165 mm	other	227-S2	227-S4

X = Injector arm
Y = Detector arm
S = Span (distance between arms)




Note
NOCI = not on-column injection
HOCl = heated on-column injection
OCI = on-column injection

Glass Columns - Silane Treated



For Shimadzu Instruments - Fittings are not supplied with Shimadzu columns

Shimadzu Instrument	Drawing (not to scale)	Column Dimensions	Length	3mm ID Item No.
GC-4BM, 14A/B, GC-6A, -7AG, -9A, 15A, 16A		X= 335 mm	1.7 m	284-1.7M3
		Y= 282 mm	2.5 m	284-2.5M3
		S= 40 mm	2.0 m	284-2.0M3
			other	284-S3
17A, 2010, 2014		X= 192 mm	1.7 m	288-1.7M3
		Y= 192 mm	2.0 m	288-2.0M3
		S= 205 mm		

For Thermo Instruments

Thermo Instruments Instrument	Drawing (not to scale)	Column Dimensions	Length	2mm ID Item No.	4mm ID Item No.
Trace 2000, NOCI Left to Left Det. Rt. to Rt. Det.		X= 180 mm	0.9 m	297-0.9M2	297-0.9M4
		Y= 200 mm	1.8 m	297-1.8M2	297-1.8M4
		S= 135 mm	2 m	297-2M2	297-2M4
			other	297-S2	297-S4
Trace 2000, NOCI Left to Rt. Det. Rt. to Left Det.		X= 200 mm	0.9 m	296-0.9M2	296-0.9M4
		Y= 175 mm	1.8 m	296-1.8M2	296-1.8M4
		S= 78 mm	2.0 m	296-2M2	296-2M4
			other	296-S2	296-S4
Trace 2000, OCI Left to Rt. Det. Rt. to Left Det.		X= 245 mm	1.8 m	298-1.8M2	298-1.8M4
		Y= 203 mm	2.0 m	298-2M2	298-2M4
		S= 78 mm			

For Varian Instruments

Varian Instrument	Drawing (not to scale)	Column Dimensions	Length	2mm ID Item No.	4mm ID Item No.
3700,4400, 4600, 6000 Vista Series, FID, ECD TSD		X= 236 mm	3ft	248-032	248-034
		Y= 203 mm	4ft	248-042	
		S= 140 mm	6ft	248-062	248-064
			2m	248-2M2	248-2M4
			8ft	248-082	
			10ft	248-102	248-104
			other	248-S2	248-S4
3300, 3400, 3600 Inj. A to Det. A		X= 203 mm	3ft	255-032	255-034
		Y= 236 mm	6ft	255-062	255-064
		S= 140 mm	other	255-S2	255-S4
Inj. B to Det. B		X= 236 mm	3ft	256-032	256-034
		Y= 203 mm	6ft	256-062	256-064
		S= 173 mm	other	256-S2	256-S4

Solid Supports and Adsorbents

Chromosorb® Supports

These are based on diatomaceous earths and are the traditional supports for gas chromatography. To alter their characteristics or make them more inert, several treatments have been developed. These include

- **NAW** (non-acid washed) - untreated
- **AW** (acid washed) - hydrochloric acid washed to remove minerals and reduce activity
- **AW-DMCS** (acid washed; Dimethyldichlorosilane treated) - silanization greatly reduces surface activity and peak tailing
- **HP** (high performance) - proprietary treatment featuring superior inertness

Chromosorb W

This is the most popular grade for GC work. It is white and relatively non adsorptive. Flux calcined. Care is needed when preparing packings and filling columns to avoid damaging the support. Maximum loading of stationary phase is 15% by weight. 100g/ bottle



Chromosorb W

Chromosorb W Type	Weight, g	Mesh Size		
		60/80	80/100	100/120
NAW	100	480001	480002	480013
AW	100	480003	480004	480005
AW-DMCS	100	480006	480008	480009
HP	100	480010	480011	480012

Chromosorb P

Chromosorb P is orange-pink and relatively hard. It is a calcined diatomite made from firebrick. Its structure makes it sturdy and capable of adsorbing more liquid phase than other supports. Maximum loading of stationary phase is 30% by weight. 100g/ bottle



Chromosorb P

Chromosorb P Type	Weight, g	Mesh Size			
		45/60	60/80	80/100	100/120
NAW	100	482001	482002	482003	
AW	100		482005	482006	
AW-DMCS	100		482008	482009	482010

Supports and Adsorbents

Chromosorb Century Series®

These are synthetic cross-linked porous polymer supports that are often used alone; they provide separations based on their pore sizes and adsorptive effects. 50g/ bottle

- **Chromosorb 101** - Proprietary polymer. For aqueous samples - alcohols, fatty acids, glycols, ethers, esters, vinyl chloride, acrylonitrile and chlorinated solvents
- **Chromosorb 102** - Styrene/ divinylbenzene. For gases and low-boiling compounds.
- **Chromosorb 103** - For amines and basic compounds
- **Chromosorb 105** - Polyaromatic resin. Separates formaldehyde from water and methanol; acetylene from lower hydrocarbons, good for most organics < 200°C
- **Chromosorb 106** - Styrene-divinylbenzene; non-polar. Gases and low-boiling compounds
- **Chromosorb 107**, **108** - Moderately polar acrylic ester resin. 107 has a surface area between 400 - 500 m²/g; 108 between 100 - 200 m²/g. 107 is preferred for trapping and analyzing vinyl acetate; 108 is for gases and aqueous compounds

Chromosorb Century Series		Mesh Size		
Type	Weight, g	60/80	80/100	100/120
101	50	483601	483701	483801
102	50	483602	483702	483802
103	50	483603	483703	483803
105	50	483605	483705	483805
106	50	483606	483706	483806
107	50	483607	483707	483807
108	50	483608	483708	483808



Chromosorb 101

Molecular Sieves

Molecular sieves are synthetic aluminosilicates of sodium, potassium or calcium, of various pore sizes. The pores are precisely uniform in size. Molecular Sieves are used to separate the fixed gases. CO, CH₄, O₂ and Ar are easily separated at room temperature. 100g/ bottle

Support	40/60	Mesh Size		
		60/80	80/100	100/120
Molecular Sieve 4A, 100 g	485320	485323	485326	485330
Molecular Sieve 5A, 100 g	485333	485336	485339	485442
Molecular Sieve 13X, 100g	485445	485448	485451	485454



Molecular Sieve 5A

Tenax®

Tenax is a porous polymer based on 2,6-diphenyl-p-phenylene oxide, used both as a column packing material and for trapping volatiles and organics. It is especially suited for high boiling liquids because of its thermal stability; the maximum temperature limit is 375°C. Tenax is an excellent trapping material since it adsorbs volatiles at room temperature, and efficiently desorbs the same materials at 300°C. Surface area 35 m²/g. Requires preconditioning.

Support	Weight, g	Mesh Size		
		35/60	60/80	80/100
Tenax	5	482021	482022	482023
Tenax	25	482025	482026	482028
Tenax	100	482029	482030	482031

Packed Column GC

Stationary Phases

Description	Sugg. Subst.	Sol	McReynold's					°C Min/Max	Item No.	Qty
			1	2	3	4	5			
Apiezon L	C		32	22	15	32	42	20/300	010017	50 g
Bentone 34	T,C							20/200	010030	50 g
Carbowax 400 (G20)	T		343	653	430			20/100	010051	*
Carbowax 1000 (G14)	M		347	607	418	626	589	40/150	010059	50 g
Carbowax 1500	C							40/200	010061	*
Carbowax 1540	C		371	639	453	666	641	40/200	010063	50 g
Carbowax 4000 (G15)	C		317	545	378	578	521	60/200	010065	50 g
Carbowax 20M (G16)	C		322	536	368	572	510	60/200	010069	50 g
Carbowax 20M Terephthalic acid (TPA)(G25)	C		321	537	367	573	520	60/250	010071	50 g
Dexsil 300 (G33)	T		47	80	103	148	96	20/450	010093	5 g
Dexsil 400	T		59	114	140	187	173	20/450	010094	5 g
Dexsil 410	T		85	165	169	242	180	20/450	010095	5 g
Dibutyl phthalate	M							-20/100	010099	*
Diethylene glycol adipate (DEGA)	A		378	603	460	665	658	20/190	010103	25 g
Diethylene glycol succinate (DEGS) (G4)	A		496	746	590	837	835	20/200	010105	25 g
Di (2-Ethylhexyl) phthalate (Octoil)	M		135	254	213	320	235	20/150	010107	50 g
Di (2-Ethylhexyl) sebacate (Octoil S) (G12)	A		72	168	108	180	125	-20/125	010109	50 g
Diisodecyl phthalate (DIDP) (G24)	A		84	173	137	218	155	-20/150	010115	50 g
Dimethylsulfolane	M							20/50	010129	10 g
Dinonyl phthalate	A		83	183	147	231	159	20/150	010133	25 g
Dioctyl phthalate	A		92	186	150	230	167	-20/100	010135	50 g
Ethylene glycol adipate	T		372	577	455	658	619	100/200	010145	50 g
FFAP	OV-351	C	340	580	397	602	627	50/250	010156	25 g
Glycerol	M							20/100	010161	50 g
Halocarbon Oil 14-25	C							20/150	010163	50 g
Igepal CO-880	C		259	461	311	482	426	100/200	010179	50 g
Kel-F Oil No. 10	A							20/100	010189	50 g
B, B-Oxydipropionitrile	M							20/100	010221	*
Phenyldiethanolamine (G12)	A							0/150	010225	*
Polyphenyl ether (5 rings) OS-124	T		176	227	224	306	283	20/200	010233	*
Polyphenyl ether (6 rings) OS-138	T		182	233	228	313	293	0/250	010235	*
Sebaconitrile	C							0/75	010244	*
Silar (see Silicone Silar)										
Silicone DC-200, 350cstk	OV-101	T	16	57	45	66	43	20/250	001010	50 g
Silicone DC-200, 12,500cstk	OV-101	T	16	57	45	66	43	0/200	001012	50 g
Silicone DC-550 (G28)	OV-101	T	81	124	124	189	145	20/225	001015	50 g
Silicone DC-710	OV-101	A	107	149	153	228	190	20/225	001018	50 g
Silicone DC-QF-1	OV-210	A	144	233	355	463	305	20/250	001021	50 g
Silicone GE-SE-30	OV-1	T	15	53	44	64	41	50/300	001023	50 g
Silicone GE-SE-52	OV-73	T	32	72	65	98	67	50/300	001025	50 g
Silicone GE-SE-54	OV-73	T	33	72	66	99	67	100/300	001026	50 g
Silicone GE-SF-96	OV-101	T	12	53	42	61	37	20/250	001027	*
Silicone GE-XE-60	OV-225	A	204	381	340	493	367	20/250	001029	*

McReynold's Code

1 = Benzene
 2 = Butanol
 3 = 2-Pentanone
 4 = Nitropropane
 5 = Pyridine

Solvent Code

A = Acetone
 T = Toluene
 C = Chloroform
 M = Methanol

*Available in packed columns only.

Packed Column GC

Stationary Phases

Description	Sugg. Subst.	Sol	McReynold's					°C Min/Max	Item No.	Qty
			1	2	3	4	5			
Silicone OV-1 (Dimethyl) (G2)		T	16	55	44	65	42	100/350	001041	10 g
Silicone OV-3 (Methyl 10% Phenyl)		A	44	86	81	124	88	20/350	001042	25 g
Silicone OV-7 (80% Methyl 20% Phenyl)		A	69	113	111	171	128	20/350	001043	25 g
Silicone OV-17 (50% Methyl 50% Phenyl)		A	119	158	162	243	202	20/350	001045	25 g
Silicone OV-25 (Phenylmethyldiphenyl) (G17)		A	178	204	208	305	280	20/350	001047	10 g
Silicone OV-61 (Diphenyldimethyl)		A	101	143	142	213	174	20/250	001048	10 g
Silicone OV-73 (Diphenyldimethyl Gum) (G27)		T	40	86	76	114	85	20/350	001049	10 g
Silicone OV-101 (Dimethyl) (G1)		T	17	57	45	67	43	20/350	001050	20 g
Silicone OV-105 (Cyanopropyldimethyl)		A	36	108	93	139	86	20/250	001051	10 g
Silicone OV-202 (Trifluoropropylmethyl) (G6)		T	146	238	358	468	310	0/275	001052	10 g
Silicone OV-210 (Trifluoropropylmethyl)		T	146	238	358	468	310	20/275	001053	25 g
Silicone OV-215 (Trifluoropropylmethyl Gum)		T	149	240	363	478	315	20/275	001057	10 g
Silicone OV-225 (Cyanopropylmethyl phenyl)		A	228	369	338	493	386	20/250	001054	10 g
Silicone OV-275 (Dicyanoallyl)		A	629	872	763	1106	849	20/275	001055	5 g
Silicone OV-330 (Carbowax-Silicone)		A	222	391	273	417	368	30/250	001056	5 g
Silicone OV-351 (Replaces FFAP)		T	335	552	382	583	540	50/250	001058	10 g
Silicone OV-1701 (Dimethylphenylcyano)		A	67	170	153	228	171	20/325	001059	3 g
Silicone Silar 5CP		T	316	494	637	531		50/275	001083	5 g
Silicone Silar 10CP		T	520	757	660	941	890	50/276	001086	5 g
Silicone UCC W-982 (G9)	OV-1	T						0/300	001035	*
SP-2100	OV-101	T	17	57	45	67	43	20/350	-	
SP-2250	OV-17	A	119	158	162	243	202	20/350	-	
SP-2401	OV-202	T	146	238	358	468	310	0/275	-	
SP-1000	OV-351	T	332	555	393	583	546	50/250	-	
SP-2310	CRS-7CPT		440	637	605	840	670	50/275	-	
SP-2330	CRS-9CPT		490	725	630	913	778	50/275	-	
SP-2300	CRS-5CPA		316	495	446	637	530	50/275	-	
SP-2340	CRS-10C	T	520	757	659	942	800	50/275	-	
Squalane		T	0	0	0	0	0	20/150	010325	*
Squalene		T	152	341	238	329	344	20/150	010327	*
THEED		M	463	942	626	801	893	20/125	010347	*
Triethanolamine		M						25/75	010355	*
1,2,3-Tris (2-Cyanoethoxy) Propane (TCEP)		M	594	857	759	1031	917	29/150	-	
Triton X-305		A	262	467	314	488	430	20/250	010365	*
UCON 50-HB-280-X		M	177	362	227	351	302	20/200	010381	50 g
UCON 50-HB-2000		A	202	394	253	392	341	20/200	010383	50 g
UCON 50-HB-5100		M	214	418	278	421	375	20/200	010385	50 g
Versamid 900								190/275	010389	*

McReynold's Code

1 = Benzene
 2 = Butanol
 3 = 2-Pentanone
 4 = Nitropropane
 5 = Pyridine

Solvent Code


A = Acetone
 T = Toluene
 C = Chloroform
 M = Methanol

*Available in packed columns only.

Packed Column GC

Solid Supports and Adsorbents

HayeSep® Porous Polymers



HayeSep D

These polymer supports are suitable for a wide range of separations. HayeSep polymers are direct replacements for the equivalent type of Porapak®; for example HayeSep Q can be directly substituted for the equivalent developed with PoraPak Q. Since there are virtually no chemically active sites in HayeSep polymers, silylation is not required.

HayeSep D - This is a high-purity divinylbenzene polymer with 80% highly-crosslinked DVB. It combines high surface area with a high operating temperature. HayeSep D polymers offer superior separation characteristics for light gases; significant separation abilities include the separation of CO and CO₂ from room air at ambient temperatures and the separation of acetylene prior to other C₂'s. It is highly recommended for the separation of water and hydrogen sulfide. HayeSep D is available in three different porosities with surface areas from 774 to 800 m²/g. This range is especially useful for certain difficult separations; for example water elutes before ethane with DB, but after ethane with Dip. Relative polarity for other HayeSep polymers is shown below; 1 is least polar and 10 is most polar.

See HayeSep Chromatograms page 40 - 41.

Type	Composition	Max. Temp*	Surface Area	Polarity**
A	DVB/EGDM	165°C	526 m ² /g	7
B	DVB/PEI	190	608	8
C	DVB/ACN	250	442	6
D	DVB, 80%	290	803	1
D _B	DVB, 80%	290	781	1
D _{IP}	DVB, 80%	290	774	1
N	DVB, EGDM	165	405	9
P	DVB, Styrene	250	165	3
Q	DVB, 60%	275	582	2
R	DVB/NVP	250	344	5
S	DVB/VP	250	583	4
T	EGDM	165	250	10

* Oxygen-free carrier gas

** 1 = lowest, 10 = highest

Polymer Materials:

- ACN - Acrylonitrile
- EGDM - Ethyleneglycoldimethacrylate
- NVP - N-Vinyl-2-Pyrrolidinone
- PEI - Polyethyleneimine
- VP - 4-Vinylpyridine
- DVB - Divinyl Benzene

HayeSep Polymers - 75 cc/ bottle

Type	Mesh Size		
	60/80	80/100	100/120
A	483021	483022	483023
B	483024	483025	483026
C	483028	483029	483030
D	483031	483032	483033
D _B	483041	483042	483043
D _{IP}	483034	483035	483036
N	483001	483002	483003
P	483004	483005	483006
Q	483014	483015	483016
R	483008	483009	483010
S	483011	483012	483013
T	483018	483019	483020

Solid Supports and Adsorbents

Porapak[®]

Porapak column packing materials consist of spherical polymer beads. These packings are chemically and physically stable with consistent particle size, porosity and surface area.

Versatility for Special Applications

Porapak packing is a versatile polymer material available in different chemistries which increase in polarity and vary in retention properties for the efficient analysis of many compounds. These packing materials offer many advantages:

- Polar molecules elute quickly making analysis simple and straightforward
- No liquid phase eliminates stationary phase bleed
- No surface hydroxyl groups eliminates tailing of highly polar molecules orof large sample loads for trace analysis and preparative GC
- Stable baselines at high temperatures



Porapak Q

Porapak Retention Data - minutes*

Compound	P	Q	R	S	N	T
Air	.23	.24	.24	.24	.24	.24
Carbon Dioxide	.36	.65	.60	.66	.90	1.26
Nitric Oxide	.27	.30	.28	.26	.28	.28
Nitrous Oxide	.42	.76	.66	.81	.94	1.22
Methane	.26	.35	.32	.38	.35	.39
Ethylene	.45	1.15	.90	1.01	1.22	1.32
Acetylene	.53	1.15	1.16	1.20	2.00	2.96
Ethane	.54	1.15	1.16	1.36	1.50	1.54
Chlorodifluoromethane	1.22	4.32	5.61	5.61	10.99	13.88
Acetic Acid	.39	1.08	2.14	3.68	4.18	5.21
Propionic Acid	.65	2.42	4.75	9.05	9.11	10.68
Methanol	.24	.35	.39	.48	.61	.69
Ethanol	.30	.60	.63	.76	1.21	1.20
Benzene	.87	2.76	2.56	3.00	4.00	3.74
Toluene	1.33	5.72	5.31	6.41	8.54	7.61
Ethyl Benzene	2.09	11.60	10.69	12.98	17.58	15.05
o-Xylene	2.52	13.48	12.53	15.32	15.32	18.05
Methylene Chloride	.41	.98	.94	1.09	1.48	1.67
Carbon Tetrachloride	.75	2.87	2.42	2.98	3.49	3.04
Methyl Acetate	.38	1.01	.97	1.11	1.74	1.75
Ethylene Glycol	.94	2.39	3.36	4.19	9.34	12.43
Glycerol	4.92	19.15	-	-	-	-

*Operating conditions: 1m x 2.3 mm ID stainless steel column, 80/100 mesh, 175°C, 25 mL/min helium flow; FID detector.

Porapak Packings - Sold in bottles of approximately 75 cc

Type	Weight, g	Mesh Size		
		50/80	80/100	100/120
P	20	489021	489022	489023
Q	26	489054	489055	489066
R	24	489088	489009	489010
S	26	489011	489012	489013
N	29	489001	489002	489003
T	31	489018	489019	489020
QS*	26	489014	489015	489016
PS*	20	489004	489005	489006

*Silanized

Packed Column GC


HayeSep® Applications

HayeSep columns are available in glass or metal, for most instruments. HayeSep porous polymers allow a range of separations still difficult to achieve with capillary columns.

HayeSep® Q For Sulphur Gas and Hydrocarbons

Column: 8' x 1/8" packed with HayeSep®, 80/100 mesh
 Column Temp: 90°C
 Flow: He 30cc/min.


1. Air
2. Methane
3. Carbon Dioxide
4. Ethylene
5. Acetylene
6. Ethane
7. Hydrogen sulfide
8. Carbonyl Sulfide
9. Propylene
10. Propane
11. Propadiene



HayeSep® B for Amines

Column: 5' x 1/8" SS packed with HayeSep® B 80/100 mesh
 Column Temp: 140°C up to 190°C at 16°C/min.
 Injector Temp: 150°C
 Detector: P.E. 900 T.C. 175 ma. 180°C
 Attenuation x 8
 Flow: He 30cc/min.
 Sample: 0.2 ml with on-column injection


1. Air
2. Water
3. Methylamine
4. Dimethylamine
5. Trimethylamine
6. Ethylene diamine



HayeSep® R for Trace Water Analysis

Column: 9' x 1/8" Ni packed with HayeSep® R, 80/100 mesh
 Column Temp: 118°C
 Flow: He at 30cc/min.
 Detector: Varain TC with Bendix On-Line analyzer
 Sample: 10 microliter, Ethyl chloride

1. Air
2. Water 12ppm
3. Hydrogen Chloride
4. Ethyl Chloride




HayeSep® C for Trace Ammonia in Water

Column: 5' x 1/8" SS HayeSep® 80/100 mesh
 Flow: He 30cc/min. 1ml sample
 Column: 115°C
 Injection: 140°C
 Detector: 180°C P.E. 900 T.C. 225 ma
 Attenuation x 2
 Chart speed: 1 cm/min.

- Trace 1
1. Air
 2. Carbon Dioxide
 3. 300 ppm Ammonia
 4. Balance Water


- Trace 2
1. Air
 2. Carbon Dioxide
 3. 100 ppm Ammonia
 4. Balance Water









HayeSep® T for Formaldehyde

Column: 5' x 1/8" packed with HayeSep® T, 100/120 mesh
 Column Temp: 132°C
 Injector Temp: 165°C
 Detector: P.E. 900 T.C. 175 ma, 180°C
 Attenuation x 32
 Flow: He 30 cc/min.
 Sample: 0.2 ml

1. Air
2. Formaldehyde
3. Water
4. Methanol



HayeSep® Applications



Packed Column GC

Custom-Coated Packings



444444

CRS provides custom-coated GC packings made to your specifications, using any of the stationary phases listed on page 36 - 37, and all common solid supports. The price of custom packings is based on the support used. Expensive stationary phases are subject to surcharge. To order, please specify the following:

- **Item No. 444444**
- **Solid Support** - ex. Chromosorb W-HP
- **Stationary phase** - ex. Carbowax 20M
- **Loading percent**- 15% maximum for Chromosorb W; 30% maximum for Chromosorb P
- **Mesh Size** - ex. 80/100
- **Quantity** - specify 20 g, 50 g, or 100 g

Support Guide

Support	Mesh Sizes Available		
Chromosorb® W-HP (higher quality W-AW-DMCS)	60/80	80/100	100/120
Chromosorb G-HP	60/80	80/100	100/120
Chromosorb 750 (newer version of W-HP)	60/80	80/100	100/120
Chromosorb 101, 102, 103, 105, 106, 107, 108	60/80	80/100	100/120
Chromosorb P (NAW, AW, HMDS, AW-DMCS)	60/80	80/100	100/120
Chromosorb G (NAW, AW, AW-DMCS)	60/80	80/100	100/120
Chromosorb W (NAW, AW, HMDS, AW-DMCS)	60/80	80/100	100/120
Chromosorb T	30/60	40/60	
HayeSep® (Type A through T)	60/80	80/100	100/120
Porapak®	50/80	80/100	

Stationary Phase Surcharge Table

The following phases are subject to surcharge:

- EGSS-X
- OV-225
- OV-275
- OV-330
- OV-1701
- Polyphenylether (6-ring)
- Silars
- Dexsils
- Tenax®

Dual phases are subject to surcharge.

Packing Accessories

Glass Wool; Quartz Wool

CRS offers five kinds of glass or quartz wool, for nearly all chromatographic needs. These are recommended for use in injection liners, and as end plugs in packed columns. Glass wool is sold in 50 g quantities; quartz wool is sold in 10 g quantities.

- **Untreated** - soft, fine-denier Pyrex® wool
- **Silane Treated** - treated with dimethyldichlorosilane (TMCS) to yield an inert material for all-around use
- **Phosphoric Acid Treated** - recommended for analysis of acidic compounds such as free acids, phenols and glycols
- **HMDS Treated** - a superior method of deactivation, using hexamethyldisilazane. Pesticide grade 10 g
- **Quartz Wool** - made from high purity quartz with low metal-oxide impurities. Ideal for use with injection port liners. Fibers have 9 um nominal diameter



Glass Wool

Description	Item No.
Untreated Glass Wool	486315
Phosphoric Treated Glass Wool	486316
Silane Treated Glass Wool	486318
HMDS Treated Glass Wool, 10 g	486319
Quartz Wool	486320

FID Cleaning Kit

Use this kit to clean both the jet tip and collector. It includes three jet reamers 0.016", 0.019" and 0.024" OD; one mini-drill bit 0.0102" OD, two wire mini-brushes, one brass and one stainless steel; and a dual ended handle for the reamers and brushes.

Description	Item No.
FID Cleaning Kit	205220



205220

Glass Wool Inserter

Used for inserting glass wool into packed columns. Simply wrap a small piece of glass wool around the tip to push it into the column.

Description	Item No.
Glass Wool Inserter	486310



486310

Glass Wool Puller

This is used for pulling glass wool out of packed columns, traps or injection port liners. It has a 1 1/4" spiral barb cut into the tip, which acts as a ratchet. For 1/4" to 2 mm IDs.

Description	Item No.
Glass Wool Puller	486313



486313

Packed Column GC

Accessories



GC End Caps

End Caps

Flexible and soft, to easily fit over the outside diameter of tubing. They are useful for sealing GC columns for storage.

Description	Qty.	Item No.
GC End Caps, 1/16"	100/pk	205121
GC End Caps, 1/8"	100/pk	205122
GC End Caps, (3/16") 5 mm	100/pk	205123
GC End Caps, 1/4"	100/pk	205124



Vibrograver

Vibrograver

Permanently engraves all hard and semi-hard surfaces; recommended for making column tags and tools. The etched information cannot be erased, wiped off, or rubbed out. Comes with a carbide tip.

Description	Voltage	Item No.
Vibrograver	110/120	486211
Vibrograver	230/240	486212
Spare Diamond Tip		486213
Spare Carbide Tip		486214




205250

Packed Column Change Kit

The Column Change Kit includes three ratchet wrenches, 7/16", 1/2" and 9/16"; one Imp® Tubing cutter, and one roll of 12.5 mm PTFE tape. The ratchet wrenches are ideal for tight spaces; the ratchet tips can be opened, then snapped closed around the fittings.

Description	Item No.
Column Change Kit	205250



486111

Column Tags, Aluminum

For packed columns; also useful for labeling computer cables, gas lines, etc. Aluminum tags are easy to inscribe and provided permanent record that cannot be erased. 22 mm x 41 mm.

Description	Qty.	Item No.
Packed Column Tags	100/pk	486111

To Order a Packed Column

To Order a Packed Column



1. Catalog Number

- Specifies the type of material - glass or stainless steel
- Specifies the length of the material - 3ft, 6ft, etc
- Specifies the shape of the column - which instrument it fits

2. Which support is required? The most common are listed below:

- Chromosorb Century Series[®]
- Chromosorb W
- Chromosorb P
- Molecular Sieve
- HayeSep[®]
- Porapak[®]

3. Specify the type of support

- Chromosorb Century Series - 101, 102, 103, 105, 106, 107, 108
- Chromosorb W - W-NAW, W-AW, W-AW-DMCS, W-HP
- Chromosorb P - P-PAW, P-AW, P-AW-DMCS
- Molecular Sieve - 4A, 5A, 13X
- HayeSep - A, B, C, D, N, P, Q, R, S, T
- Porapak - P, Q, R, S, N, T

4. What mesh size is required?

- Chromosorb - 60/80, 80/100, 100/120
- Molecular Sieve - 40/60, 60/80, 80/100, 100/120
- HayeSep - 60/80, 80/100, 100/120
- Porapak - 50/80, 80/100, 100/120

5. What type of Stationary Phase is needed (if any)? Up to 2 phases

- Chromosorb Century Series, Molecular Sieve, HayeSep, Porapaks do not require a stationary phase
- Chromosorb W and Chromosorb P require a stationary phase

6. What % loading is needed on the stationary phase?

- Chromosorb W - maximum loading is 15%
- Chromosorb P - maximum loading is 30%

Surcharge Table: Surcharges apply to the following expensive phases

Stationary Phase Surcharges

- Dexsils
- EGSS-X
- OV-225
- OV-275
- OV-330
- OV-1701
- Polyphenylether
- Silars

Support Surcharges

- Chromosorb T
- Carbon Molecular Sieves
- HayeSep
- Porapak
- Chromosorb 101 - 109
- Tenax[®]
- Carbopacks