

# ULTRAPURE WATER SYSTEMS

## Integrity series

Crystal Integrity systems produce ultrapure water directly from tap water. Ultrapure water produced by the Crystal Integrity system complies with the requirements for Grade 1 water in accordance with relevant standards (ISO 3696, ASTM, CLSI).

The Crystal Integrity systems are recommended for laboratories with relatively low consumption of ultrapure water (about 5 litres per day).

Crystal Integrity systems are available in the following configurations:

- **Integrity HPLC System** (P/N CB-1503) produces water with very low organic carbon (TOC) content. This meets the requirements for liquid chromatography methods. Integrity HPLC water can also be used for some microbiological and molecular biology applications.

- **Integrity Bio System** (P/N CB-1505) produces water with very low organic and RNase/DNase content. This is intended for use in molecular biology, including RNase-sensitive applications.

Crystal Integrity water is intended for the most demanding applications including, but not limited to:

- inorganic trace analysis
- liquid chromatography
- cell culture
- molecular biology

Ultrapure water produced by an Integrity system, with resistivity of 18.2 MegaOhm\*cm (0.055 µS/cm), exceeds requirements of all relevant standards (ISO 3696 Grade 1, ASTM Type I, CLSI Type I).

Purified water is collected in an integrated storage tank. The recirculation system ensures consistent quality of water and reduces total organic carbon (TOC) to very low levels: <2ppb.

Space-saving, integrated, 5L tank is made of high density polyethylene to minimize contamination. The storage tank is accessible for cleaning.

All Integrity systems have a controller with a graphic LCD display for water quality indication. The LCD display provides all necessary information about system status, remaining pre-filter life, and polishing module performance. The smart polishing module monitoring system provides a reduction in running costs. Users are instructed to replace the polishing module only when the module is close to the end of its service life.

If a Bluetooth option (P/N 10103) is installed, the controller can be detached. The Integrity system can be controlled remotely for more convenient operation.

All cartridges and filters are easily accessible and no tools are required to replace them. The Integrity system can be installed on a laboratory bench or mounted on the wall to save space.

The Crystal Integrity has important safety functions, including:

- tank filling control
- tap water pressure control

## Applications

Application	Integrity HPLC	Integrity Bio
Reagent preparation	●	●
Ion chromatography	●	●
ICP-MS	●	●
Atomic absorption spectrophotometry with graphite atomizer	●	●
ICP-OES	●	●
HPLC	●	●
Gas chromatography	●	●
Total Organic Carbon measurement	●	●
Flow cytometry	–	●
Cell and tissue culture	–	●
Molecular biology	–	●



## Components

Part number	Accessory	Integrity HPLC	Integrity Bio
10019	Pre-filter set for chlorine and organics reduction	+	+
	Reverse osmosis module	+	+
10029	Polishing module	+	+
10105	Photo-oxidation module	+	+
10012	Point-of-use microfilter	+	–
10109	Point-of-use ultrafiltration module	–	+
10102	UV sterilization module	–	+
10106	Integrated TOC monitor	option	–
10103	Removable controller with Bluetooth module	option	option

## Consumables

Part number	Description	Replacement criteria	Comments
10019	Replacement pre-filter set, Integrity	Filter life counter is zero or the filter is clogged	
10029	Replacement polishing module	„DI Err“ message is shown, or water conductivity while dispensing is consistently > 0.1 µS/cm	
10011	Replacement sterilization UV bulb	As required (on average – every 3 years)	„Bio“ systems only
10018	Replacement photooxidation UV bulb	As required (on average – every 3 years)	„HPLC“ and „Bio“ systems
10012	Replacement 0.22 µm dispense microfilter	Every 6–12 months	„Trace“ and „HPLC“ systems
10120	Replacement ultrafilter	Every 6–12 months	„Bio“ systems only

## Flow diagrams



## Specifications

Purified water specifications	Integrity HPLC	Integrity Bio
Ultrapure (Grade 1) water resistivity	18.2 MΩ x cm	18.2 MΩ x cm
Ultrapure (Grade 1) water conductivity	0.055 µS/cm	0.055 µS/cm
TOC	<2 ppb	<2 ppb
RNase	–	<0.01 ng/mL
DNase	–	<4 pg/µL
Bacteria	< 1 cfu/mL	< 1 cfu/mL
Endotoxins	<0.15 EU/mL	< 0.001 EU/mL
Particles > 0.22 µm	<1/ per mL	<1/ per mL
Nominal flow, pure water (to storage tank)	4 L/h	4 L/h
Nominal dispense flow, ultrapure water	2 L/min	2 L/min
Polishing module life	1 m <sup>3</sup>	1 m <sup>3</sup>
Recovery	>30 %	>30 %
Dimensions (WxDxH), cm	50x40x60	50x40x60
Integrated storage tank	5 L	5 L
Feed water pressure	1.5 – 5 bar	1.5 – 5 bar
Feed water conductivity	< 900 µS/cm	< 900 µS/cm



## Ordering Information

Model	Part number
Crystal Integrity HPLC	CB-1503
Crystal Integrity Bio	CB-1505