

GsBP-PONA

Excellent column for gasoline, natural gas, petroleum gases, sulfurs, amines, alcohols and other oxygenates

- 100% Poly(dimethylsiloxane), crosslinked and bonded
- Solvent rinseable, high temperature limit
- Low column bleed, inert and neutral surface, very high plate numbers
- 50m and 100m versions
- Non-polar

Similar Phases: HP-PONA, DB-Petro, Rtx-PONA

GsBP-PONA Ordering Info :

ID(mm)	Length(m)	Film(um)	Temperature°C	P/N
0.2	50	0.50	-60 to 325/350	9002-PONA
0.25	100	0.50	-60 to 325/350	9006-PONA