

DURACAL pH buffers

Simple handling for professional results:



Step 1 Open bottle



Step 2 Fill calibration compartment



Step 3 Calibrate electrode



Step 4 Empty calibration compartment

pH buffers

pH Value	Accuracy	Stability (in months)	Certified by	Packaging unit	PN
1.09	±0.02	60	HAMILTON	500 mL	238271
1.68	±0.02	60	HAMILTON	500 mL	238272
2.00	±0.02	60	HAMILTON	500 mL	238273
3.06	±0.02	60	HAMILTON	500 mL	238274
4.01	±0.01 / ±0.02	24 / 60	DKD	250 mL	238317
4.01	±0.01 / ±0.02	24 / 60	DKD	500 mL	238217
4.01	±0.01 / ±0.02	24 / 60	DKD	3 x 500 mL	238917
4.01	±0.01 / ±0.02	24 / 60	DKD	5 L	238332
4.01	±0.01 / ±0.02	24 / 60	DKD	10 L	238194
4.01	±0.01 / ±0.02	24 / 60	DKD	1000 L	238895
5.00	±0.02	60	HAMILTON	500 mL	238275
6.00	±0.02	60	HAMILTON	500 mL	238276
7.00	±0.01 / ±0.02	24 / 60	DKD	250 mL	238318
7.00	±0.01 / ±0.02	24 / 60	DKD	500 mL	238218
7.00	±0.01 / ±0.02	24 / 60	DKD	3 x 500 mL	238918
7.00	±0.01 / ±0.02	24 / 60	DKD	5 L	238333
7.00	±0.01 / ±0.02	24 / 60	DKD	10 L	238188
7.00	±0.01 / ±0.02	24 / 60	DKD	1000 L	238896
8.00	±0.02	60	HAMILTON	500 mL	238277
9.21	±0.02	60	DKD	250 mL	238319
9.21	±0.02	60	DKD	500 mL	238219
9.21	±0.02	60	DKD	3 x 500 mL	238919
9.21	±0.02	60	DKD	10 L	238216
9.21	±0.02	60	DKD	1000 L	238897
10.01	±0.02	60	DKD	250 mL	238321
10.01	±0.02	60	DKD	500 mL	238223
10.01	±0.02	60	DKD	3 x 500 mL	238923
10.01	±0.02	60	DKD	10 L	238187
10.01	±0.02	60	DKD	1000 L	238898
11.00	±0.02	24	HAMILTON	500 mL	238278
12.00	±0.02	24	HAMILTON	500 mL	238279
4.01/7.00/9.21	±0.01 / ±0.02	24 / 60	DKD	500 mL, mixed	238922
4.01/7.00/10.01	±0.01 / ±0.02	24 / 60	DKD	500 mL, mixed	238924

Redox buffers

Redox value	Accuracy	Stability (monthly)	Certified by	Packing unit	PN
271 mV	±5 mV	24	None	500 mL	238228
475 mV	±5 mV	24	None	250 mL	238322
475 mV	±5 mV	24	None	500 mL	238227