

Air Generator

The *Peerless* Air Generator produces a continuous flow of purified air from the built-in air compressor. The gas

is filtered by two Activated Carbon traps. This gas is good for use as a fuel gas in GC/FID.



Features

- Built-in Air Compressor
- Noise Level <35 db
- Compact, Light Weight Design
- CE Certification

Benefits

- Replaces Air Gas Cylinder
- Convenience...Gas on demand
- Consistent Purity
- Safe...Very low volume stored at a max 5kgs/cm²

Models	Flowrate cc/min	Pressure	Voltage Vac	Electrical Consumption	Dimension (mm)			Weight kg
					A	B	C	
SPB-20	2000	5 kg/cm ²	220	175W	420	250	350	20
SPB-50	5000	5 kg/cm ²	220	275W	420	250	350	20