

PEERLESS LC (DEACTIVATED FOR BASES)

PHYSICAL CHARACTERISTICS

Phase	Silica Type	Family	% C	Pore Size Å	Particle Sizes (µm)	End-Capped Yes/No	Surface Area m ² /g	USP Code
<i>Peerless LC</i>	Type A (Deactivated for Bases)	C18	11%	120	3,5	Y	180	L-1
		C8	7%	120	3,5	Y	180	L-7
		Phenyl	5%	120	5	Y	180	L-11
		Amino	3%	100	3,5	N	280	L-8
		Diol	3%	100	5	N	280	L-20

Part No.	Range
2546521	25cm x 4.6mm x 5µm Peerless LC C18
2540521	25cm x 4.0mm x 5µm Peerless LC C18
1546521	15cm x 4.6mm x 5µm Peerless LC C18
1540521	15cm x 4.0mm x 5µm Peerless LC C18
1046521	10cm x 4.6mm x 5µm Peerless LC C18
1040521	10cm x 4.0mm x 5µm Peerless LC C18
0546521	5cm x 4.6mm x 5µm Peerless LC C18
0540521	5cm x 4.0mm x 5µm Peerless LC C18
1546321	15cm x 4.6mm x 3µm Peerless LC C18
1540321	15cm x 4.0mm ID x 3µm Peerless LC C18
1046321	10cm x 4.6mm x 3µm Peerless LC C18
1040321	10cm x 4.0mm ID x 3µm Peerless LC C18
0546321	5cm x 4.6mm x 3µm Peerless LC C18
0540321	5cm x 4.0mm x 3µm Peerless LC C18

2546522	25cm x 4.6mm x 5µm Peerless LC C8
1546522	15cm x 4.6mm x 5µm Peerless LC C8
1046522	10cm x 4.6mm x 5µm Peerless LC C8
0546522	5cm x 4.6mm x 5µm Peerless LC C8
1546322	15cm x 4.6mm x 3µm Peerless LC C8
1046322	10cm x 4.6mm x 3µm Peerless LC C8
0546322	5cm x 4.6mm x 3µm Peerless LC C8

2546525	25cm x 4.6mm x 5µm Peerless LC Phenyl
1546525	15cm x 4.6mm x 5µm Peerless LC Phenyl
1046525	10cm x 4.6mm x 5µm Peerless LC Phenyl
0546525	5cm x 4.6mm x 5µm Peerless LC Phenyl
1546325	15cm x 4.6mm x 3µm Peerless LC Phenyl
1046325	10cm x 4.6mm x 3µm Peerless LC Phenyl
0546325	5cm x 4.6mm x 3µm Peerless LC Phenyl

2546526	25cm x 4.6mm x 5µm Peerless LC Amino
1546526	15cm x 4.6mm x 5µm Peerless LC Amino
1046526	10cm x 4.6mm x 5µm Peerless LC Amino
0546526	5cm x 4.6mm x 5µm Peerless LC Amino
1546326	15cm x 4.6mm x 3µm Peerless LC Amino
1046326	10cm x 4.6mm x 3µm Peerless LC Amino
0546326	5cm x 4.6mm x 3µm Peerless LC Amino

2546527	25cm x 4.6mm x 5µm Peerless LC Diol
1546527	15cm x 4.6mm x 5µm Peerless LC Diol
1046527	10cm x 4.6mm x 5µm Peerless LC Diol
0546527	5cm x 4.6mm x 5µm Peerless LC Diol
1546327	15cm x 4.6mm x 3µm Peerless LC Diol
1046327	10cm x 4.6mm x 3µm Peerless LC Diol
0546327	5cm x 4.6mm x 3µm Peerless LC Diol

These columns are also offered in other geometries. Please enquire with us for unlisted part numbers.