

Electronic Pipettes

SoftGrip™
Pipettes

SofTouch™
Electronic Pipettes





AdvanTip™
Pipette Tips

SoftAide™
Pipette Controller





SoftAide™
Quik Dispenser

Replacement Parts
and Accessories




SofTouch™ Single Channel Electronic Pipettes, 0.2 µL - 5 mL

Part No.	Volume	Increments (µL)	Color		At 10% of Pipette Volume		At 100% of Pipette Volume	
					Accuracy within ±%	Precision within +%	Accuracy within ±%	Precision within +%
8707-01XX	0.2-10 µL	0.1	Purple		2.50	1.50	0.90	0.50
8707-02XX	5-100 µL	1.0	Violet		2.00	1.00	0.40	0.15
8707-03XX	50-1000 µL	5.0	Sky Blue		1.50	1.50	0.40	0.15
8707-04XX	100-5000 µL	50.0	Graphite		0.80	0.30	0.50	0.15

SofTouch Multi-Channel Electronic Pipettes, 8-Channel, 0.2 µL - 1200 µL

Part No.	Volume	Increments (µL)	Color		At 10% of Pipette Volume		At 100% of Pipette Volume	
					Accuracy within ±%	Precision within +%	Accuracy within ±%	Precision within +%
8707-05XX	0.2-10 µL	0.1	Purple		4.00	4.00	0.90	0.50
8707-06XX	5-100 µL	1.0	Violet		2.50	1.50	0.50	0.20
8707-07XX	25-250 µL	5.0	Burnt Orange		1.50	1.00	0.40	0.10
8707-08XX	50-1200 µL	10.0	Teal		4.00	0.80	0.50	0.15

SofTouch Multi-Channel Electronic Pipettes, 12-Channel, 0.2 µL - 1200 µL

Part No.	Volume	Increments (µL)	Color		At 10% of Pipette Volume		At 100% of Pipette Volume	
					Accuracy within ±%	Precision within +%	Accuracy within ±%	Precision within +%
8707-09XX	0.2-10 µL	0.1	Purple		4.00	4.00	0.90	0.50
8707-11XX	25-250 µL	5.0	Burnt Orange		1.50	1.00	0.40	0.10
8707-12XX	50-1200 µL	10.0	Teal		4.00	0.80	0.80	0.15

Replace 'XX' for regional adapter: US = U.S.A., EU = Europe, UK = United Kingdom, and JP = Japan.

Direct-connect AC/DC recharging adapter and environmentally friendly NiMH rechargeable battery included with all SofTouch Pipettes.

Optional Pipette Safe™ Tip Cone Filters to prevent sample cross-contamination are available for all SofTouch Pipettes.

Selecting the Right Hamilton AdvanTip™ Pipette Tip

Pipette Part No.	Pipette Volume	Hamilton Tip	Page No.
8707-01XX	0.2-10 µL	10 µL	21 and 22
8707-02XX	5-100 µL	100 and 200 µL	21 and 22
8707-03XX	50-1000 µL	1000 µL	21 and 22
8707-04XX	100-5000 µL	5000 µL	21 and 22
8707-05XX	0.2-10 µL x 8	10 µL	21 and 22
8707-06XX	5-100 µL x 8	100 and 200 µL	21 and 22

Pipette Part No.	Pipette Volume	Hamilton Tip	Page No.
8707-07XX	25-250 µL x 8	300 µL	21 and 22
8707-08XX	50-1200 µL x 8	1200 µL	21 and 22
8707-09XX	0.2-10 µL x 12	10 µL	21 and 22
8707-11XX	25-250 µL x 12	300 µL	21 and 22
8707-12XX	50-1200 µL x 12	1200 µL	21 and 22



SofTouch 5-Position Carousel Pipette Stand

Safe, convenient storage solution for up to five single channel SofTouch pipettes.

Part No.	Description
8706-01	Holds Five SofTouch Pipettes