

Autosampler Vials

- Borosilicate Glass- USP Type I, Class A in Clear and Amber
- Snap-Top Design
- Crimp-Top Design
- Screw-Top Design
- Shell Design
- Vials for Headspace Sampling
- Low Volume Glass Inserts
- Plastic (Polypropylene) Vials also on offer
- ISO 9002 Certified Manufacturing

Vial Accessories

- Crimpers
- Decrimpers
- Vial Closures-Assorted Caps and Seals
- Vial Racks

Properties	Clear	Amber
Linear Coefficient of Expansion	33	51
Maximum Temperature (° C)	515	535
USP Class Type	Type 1	Type 1
Light Protection	No	Yes



Screw Top Vials

Part No.

GCS-1232-08-0

GAS-1232-08-0

Standard Mouth (8 x 425) 12mm x 32mm

2ml Clear Screw Top Glass Vial (pk/100)

2ml Amber Screw Top Glass Vial (pk/100)

Closures

PKSC-08-0

Black Screw Caps (pk/100)

PKSC-08-2

Black Screw Caps with Red PTFE/Silicone Septa (pk/100)

SCS-08-2

8mm Red PTFE/Silicone Septa (pk/100)

SCS-08-3B

8mm Blue PTFE/WhiteSilicone pre-slit Septa (pk/100)



GCS-1232-09-0

GCS-1232-09-M

GAS-1232-09-M

GAS-1232-09-0

PPS-1232-09-LV5

Wide Mouth (9 x 425) 12mm x 32mm

2ml Clear Screw Top Glass Vial (pk/100)

2ml Clear Screw Top Glass Vial with write-on spot (pk/100)

2ml Amber Screw Top Glass Vial with write-on spot (pk/100)

2ml Amber Screw Top Glass Vial (pk/100)

Low Volume 500 ul capacity Polypropylene Vial (pk/100)

Closures

PBSC-09-1

Blue Screw Cap with PTFE/Red Rubber Septa (pk/100)

PBSC-09-2

Blue Screw Cap with PTFE/Silicone Septa (pk/100)

PBSC-09-3

Blue Cap with Slit PTFE/Silicone Septa (pk/100)

PBSC-09-0

Blue Open Top Screw Caps (pk/100)

PKBSC-09-3

Knurled Blue Screw Cap with Slit PTFE/Silicone Septa (pk/100)

SCS-09-3

9mm Preslit Red PTFE/White Silicone Septa (pk/100)

SCS-09-1

PTFE/Red Rubber Septa (9mm) (pk/100)

SCS-09-2

Red PTFE/White Silicone Septa (9mm) (pk/100)

309044

Blue Screw Cap with Bonded PTFE/Silicone Slit Septa (pk/100)



Standard Mouth (13 x 425) 15mm x 45mm	
GCS-1545-13-0	4ml Clear Screw Top Vial (pk/100)
GCA-1545-13-0	4ml Amber Screw Top Vial (pk/100)

Closures

SCS-13-2	13mm Red PTFE/Silicone Septa (pk/100)
PKSC-13-2	Black Cap with PTFE/Silicone Septa (pk/100)
PKSC-13-0	13mm Black Screw Cap (pk/100)

**Crimp Top Glass Vials**

Wide Mouth 12mm x 32mm	
GCC-1232-WM-0	2ml Clear Glass Crimp Top Vial (pk/100)
GCC-1232-WM-M	2ml Clear Glass Crimp Top Vial with write-on spot (pk/100)
GAC-1232-WM-0	2ml Amber Glass Crimp Top Vial (pk/100)
GAC-1232-WM-M	2ml Amber Glass Crimp Top Vial with write-on spot (pk/100)

Closures

ASCC-11-1	11mm Aluminium Cap with PTFE/Red Rubber Septa (pk/100)
ASCC-11-2	11mm Aluminium Cap with PTFE/Silicone Septa (pk/100)
ASCC-08-1	8mm Aluminium Cap with PTFE/Red Rubber Septa (pk/100)

**Headspace Vials**

GCC-2375-SM-F	20ml (23mm x 75mm) Flat Bottom Vials (pk/100)
GCC-2375-SM-R	20ml (23mm x 75mm) Round Bottom Vials (pk/100)
GCC-2346-SM-F	10ml (23mm x 46mm) Flat Bottom Vials (pk/100)
GCC-3060-SM-F	27ml (30mm x 60mm) Vial for Shimadzu (pk/100)

Closures

ASCC-20-2	20mm Tan PTFE/White Silicone Septa with Aluminium Crimp Cap (pk/100)
SCC-20-2	20mm Tan PTFE/White Silicone Septa (pk/100)
ASCC-20-2W	20mm White PTFE/White Silicone Septa with Aluminium Crimp Cap (pk/100)
SCS-18-2	18mm Blue PTFE/White Silicone Septa for Metal Screw Cap (pk/100)

**Snap-Top Vial**

Wide Mouth 12mm x 32mm	
GCN-1232-WM-0	2ml Clear Glass Snap Top Vial (pk/100)
GCN-1232-WM-M	2ml Clear Glass Snap Top Vial with write-on spot (pk/100)
GAN-1232-WM-M	2ml Amber Glass Snap Top Vial with write-on spot (pk/100)

Closures

PCNC-12-4	Clear Polyethylene Snap Caps with Solid Polyethylene Membrane Septa (pk/100)
PCNC-12-2	Clear Polyethylene Snap Caps with Red PTFE/Silicone Septa (pk/100)
PCNC-12-1	Clear Polyethylene Snap Caps with PTFE/Red Rubber Septa (pk/100)
PCNC-12-5	Clear Polyethylene Snap Cap with Moulded Star Bust Septa (pk/100)



Shell Vial

GCH-0840-08-0	8mm X 40mm Shell Vial for Waters WISP 96 Position Autosampler
GAH-0840-08-0	1.2ml Clear Glass Shell Vial without Plug (pk/100)
PCH-0840-08-0	1.2ml Amber Glass Shell Vial without Plug (pk/100)
	1.2ml Polypropylene Shell Vial without Plug (pk/100)

12mm X 32mm Shell Vials

GCH-1232-12-0	2ml Clear Glass Shell Vial without Plug (pk/100)
GAH-1232-12-0	2ml Amber Glass Shell Vial without Plug (pk/100)
PCH-1232-12-0	2ml Polypropylene Shell Vial without Plug (pk/100)

Closures-Polyethylene Snap Plugs

PCNP-08-5	8mm Clear Snap Plug with Star Bust (pk/100)
PCNP-12-5	12mm Clear Snap Plug with Star Bust (pk/100)



Inserts

Low Volume inserts for 9mm thread Vials

GCI-0630-09-4	50µl Glass Insert with bottom spring (pk/100)
GCI-0630-09-3	100µl Glass Insert with spring (pk/100)
GCI-0630-09-2	250µl Glass Insert without spring (pk/100)
GCI-0630-09-1	250µl Glass Insert with bottom spring (pk/100)
GCI-0630-09-0	350µl Flat bottom insert (pk/100)

Low Volume inserts for 8mm thread standard mouth Vials

GCI-0530-08-3	50µl Glass Insert with bottom spring (pk/100)
GCI-0530-08-2	100µl Glass Insert with bottom spring (pk/100)
GCI-0530-08-1	200µl Glass Insert without bottom spring (pk/100)
GCI-0530-08-0	250µl Flat bottom insert (pk/100)



Vial Racks

PBVR-U	Universal Rack
PBVR-08-25	25 Position Vial Rack for 8mm Vials
PBVR-12-25	25 Position Vial Rack for 12mm Vials
PBVR-15-25	25 Position Vial Rack for 15mm Vials
PBVR-08-96	96 Position Vial Snap Rack for 8mm Vial
PBVR-12-50	50 position Vial Snap Rack for 12mm Vial



Crimper & Decrimper

Electronic Crimper

- Size 11mm
- Size 20mm

Manual Crimper

- Size 8mm
- Size 11mm
- Size 20mm

Manual Decrimper

- Size 8mm
- Size 11mm
- Size 20mm

