

GsBP-5 Ordering Info

ID (mm)	Length (m)	Film (μm)	Temperature Limit ($^{\circ}\text{C}$)	P/N	
0.20	12	0.33	-60 to 325/350	0520-1203	
	17	0.33	-60 to 325/350	0520-1703	
	25	0.11	-60 to 325/350	0520-2501	
	25	0.33	-60 to 325/350	0520-2503	
	25	0.50	-60 to 325/350	0520-2505	
	50	0.11	-60 to 325/350	0520-5001	
	50	0.33	-60 to 325/350	0520-5003	
	50	0.50	-60 to 325/350	0520-5005	
	0.25	15	0.10	-60 to 325/350	0525-1501
15		0.25	-60 to 325/350	0525-1502	
15		0.50	-60 to 325/350	0525-1505	
15		1.0	-60 to 325/350	0525-1510	
30		0.10	-60 to 325/350	0525-3001	
30		0.25	-60 to 325/350	0525-3002	
30		0.50	-60 to 325/350	0525-3005	
30		1.00	-60 to 325/350	0525-3010	
60		0.10	-60 to 325/350	0525-6001	
60		0.25	-60 to 325/350	0525-6002	
60		0.50	-60 to 325/350	0525-6005	
60		1.00	-60 to 325/350	0525-6010	
0.32	15	0.10	-60 to 325/350	0532-1501	
	15	0.25	-60 to 325/350	0532-1502	
	15	0.50	-60 to 325/350	0532-1505	
	15	1.00	-60 to 325/350	0532-1510	
	25	0.17	-60 to 325/350	0532-2501	
	25	0.52	-60 to 325/350	0532-2505	
	25	1.05	-60 to 325/350	0532-2510	
	30	0.10	-60 to 325/350	0532-3001	
	30	0.25	-60 to 325/350	0532-3002	
	30	0.50	-60 to 325/350	0532-3005	
	30	1.00	-60 to 300/320	0532-3010	
	50	0.17	-60 to 325/350	0532-5001	
	50	0.25	-60 to 325/350	0532-5002	
	50	0.52	-60 to 325/350	0532-5005	
	50	1.05	-60 to 300/320	0532-5010	
	60	0.10	-60 to 325/350	0532-6001	
	60	0.25	-60 to 325/350	0532-6002	
	60	1.00	-60 to 325/350	0532-6010	
	0.53	10	2.65	-60 to 260/280	0553-1026
		10	5.00	-60 to 260/280	0553-1050
		15	1.50	-60 to 300/320	0553-1515
15		3.00	-60 to 260/280	0553-1530	
15		5.00	-60 to 260/280	0553-1550	
30		0.25	-60 to 300/320	0553-3002	
30		0.88	-60 to 300/320	0553-3008	
30		1.00	-60 to 300/320	0553-3010	
30		1.50	-60 to 300/320	0553-3015	
30		2.65	-60 to 260/280	0553-3026	
30		3.00	-60 to 260/280	0553-3030	
30		5.00	-60 to 260/280	0553-3050	
60		3.00	-60 to 260/280	0553-6030	
60		5.00	-60 to 260/280	0553-6050	