

Acidic/Neutral Drugs	Endocrine Disruptors: Phthalate Esters	Formaldehyde, 50ppb	Polybrominated Diphenyl Esters (PBDE)
Chlorinated Pesticides, EPA Method 508	EPA Method 8061 (Phthalate Ester)	Organochlorine Pesticides	Sulfur in air
CLP Standard, Semivolatile Organics	EPA method 8270	Organochlorine Pesticides II EPA Method 8081A	Trace Active Amines, 10 ng on-column
Diesel Fuel	EPA-625 Phenols	Organohalide Pesticides in water, EPA 505	Tricyclic Antipsychotics
Drug of abuse	Flavor Mixture	Organophosphorus Pesticides I EPA 8141A	Urine Drug Screen
Endocrine Disruptors: Alkyl Phenols	Food Packaging Volatiles	Phenols, I and II	US EOA Method 8270D mix

GsBP-5MS Ordering Info

ID (mm)	Length (m)	Film (μm)	Temperature Limit (°C)	P/N
0.20	12	0.33	-60 to 325/350	1520-1203
	17	0.33	-60 to 325/350	1520-1703
	25	0.33	-60 to 325/350	1520-2503
	50	0.33	-60 to 325/350	1520-5003
0.25	15	0.10	-60 to 325/350	1525-1501
	15	0.25	-60 to 325/350	1525-1502
	15	0.50	-60 to 325/350	1525-1505
	15	1.00	-60 to 325/350	1525-1510
	30	0.10	-60 to 325/350	1525-3001
	30	0.25	-60 to 325/350	1525-3002
	30	0.50	-60 to 325/350	1525-3005
	30	1.00	-60 to 325/350	1525-3010
	60	0.10	-60 to 325/350	1525-6001
	60	0.25	-60 to 325/350	1525-6002
	60	0.50	-60 to 325/350	1525-6005
	60	1.00	-60 to 320/340	1525-6010
0.32	15	0.10	-60 to 325/350	1532-1501
	15	0.25	-60 to 325/350	1532-1502
	15	0.50	-60 to 325/350	1532-1505
	15	1.00	-60 to 325/350	1532-1510
	25	0.52	-60 to 325/350	1532-2505
	30	0.10	-60 to 325/350	1532-3001
	30	0.25	-60 to 325/350	1532-3002
	30	0.50	-60 to 325/350	1532-3005
	30	1.00	-60 to 325/350	1532-3010
	60	0.10	-60 to 325/350	1532-6001
	60	0.25	-60 to 325/350	1532-6002
	60	0.50	-60 to 325/350	1532-6005
	60	1.00	-60 to 300/325	1532-6010
	0.53	0.50	-60 to 300/320	1553-3005
		0.88	-60 to 300/320	1553-3008
		1.50	-60 to 300/320	1553-3015

GsBP-5MS with Integrated Guard Column

In addition to all features of GsBP-5MS columns, the integrated guard columns provide:

- Extended column life
- Improved sample applications
- Same inertness as regular 5-MS columns
- Improved peak shapes for low boiling point compounds with presence of high boiling point compounds

Similar Phases: DB-5MS, HP-5MS, Ultra-5, SLB-5, Rxi-5MS, Rxi-5sil MS, Rxi-5HT, VF-5MS, CP-Sil 8CB low bleed MS, BPX-5, ZB-5MS, ZB-5HT, AT-5MS, Optima-5MS