

Marathon™ Septa

- Excellent Performance with Long Injection Life
- Pre-pierced for greater reliability
- Premium septum for autosamplers
- CenterGuide™ Marathon are pre-pierced for extremely reliable performance
- Maximum temperature setting 400°C

The CenterGuide versions of Marathon are pre-pierced to channel the needle through the same spot with each injection, which helps eliminate the common type of coring that occurs on the first injection. The Marathon septum is also ideal for standard manual-injection GC and GC/MS. Septa are 3 mm thick, except Shimadzu-style plug septa. Septa are packaged in glass jars.



Marathon 11mm Septa

Description	25/pk	50/pk	100/pk
5 mm Marathon CenterGuide	239587	239597	
9 mm Marathon CenterGuide	239778	239779	
7/16" (11 mm) Marathon CenterGuide	239287	239297	239888
11.5 mm Marathon CenterGuide	239787	239798	
1/2" (12.7 mm) Marathon CenterGuide	239387	239397	
17 mm Marathon	239688	239698	
Shimadzu	239488	239498	