

SoftGrip™ Pipettes – Single Channel

- **Wide Selection of Ergonomic Pipettes**
- **Completely Autoclavable**
- **Guaranteed Accuracy and Precision**

SoftGrip pipettes routinely receive the highest user satisfaction score of 10 on 11 different aspects of a pipette, including plunger and eject forces, weight, diameter, and texture. Choose from six adjustable volume and ten fixed volume SoftGrip pipettes covering 0.2 μ L to 1 mL. All pipettes have color-coded plunger buttons and labels for easy volume identification.

SoftGrip air-displacement pipettes are built from high quality, state-of-the-art materials designed to provide years of reliable service. All pipettes are completely autoclavable for easy decontamination and sterilization without the need for disassembly and recalibration.

The isolated volume ring prevents accidental volume changes. A magnified, easy-to-read micrometer ensures reproducible, accurate volume setting. Each pipette meets or exceeds industry standards for accuracy and precision when used with Hamilton Universal AdvanTip™ Precision Pipette Tips.



Manual Single Channel Pipettes

SoftGrip™ Adjustable Volume Pipettes, 0.2 µL - 1 mL

Part No.	Volume	Increments (µL)	Color	At 10% of Pipette Volume		At 100% of Pipette Volume	
				Accuracy within ±%	Precision within +%	Accuracy within ±%	Precision within +%
55019-30	0.2-2 µL	0.002	Aqua	8.00	4.00	1.20	0.60
55019-32	1-10 µL	0.02	Purple	2.50	1.20	0.80	0.40
55019-34	2.5-25 µL	0.02	Forest Green	4.50	1.50	0.80	0.20
55019-36	10-100 µL	0.2	Violet	1.80	0.70	0.60	0.15
55019-38	30-300 µL	0.2	Brick Red	1.20	0.40	0.40	0.15
55019-40	100 µL-1 mL	2.0	Sky Blue	1.60	0.50	0.30	0.12

SoftGrip Fixed Volume Pipettes, 5 µL - 1 mL

Part No.	Volume	Color	At 100% of Pipette Volume	
			Accuracy within ±%	Precision within +%
55019-01	5 µL	Steel Blue	1.40	0.75
55019-03	10 µL	Purple	0.80	0.40
55019-05	25 µL	Forest Green	0.80	0.30
55019-07	50 µL	Sandstone	0.60	0.30
55019-09	100 µL	Violet	0.50	0.20
55019-11	200 µL	Mustard	0.40	0.18
55019-13	250 µL	Burnt Orange	0.40	0.18
55019-15	300 µL	Brick Red	0.40	0.18
55019-17	500 µL	Olive Green	0.40	0.18
55019-19	1 mL	Sky Blue	0.30	0.12

Selecting the Right Hamilton AdvanTip™ Pipette Tip

Pipette Part No.	Pipette Volume	Hamilton Tip	Page No.	Pipette Part No.	Pipette Volume	Hamilton Tip	Page No.
55019-30	0.2-2 µL	10 µL	21 and 22	55019-01	5 µL	10 µL	21 and 22
55019-32	1-10 µL	10 µL	21 and 22	55019-03	10 µL	10 µL	21 and 22
55019-34	2.5-25 µL	200 µL	21 and 22	55019-05	25 µL	200 µL	21 and 22
55019-36	10-100 µL	200 µL	21 and 22	55019-07	50 µL	200 µL	21 and 22
55019-38	30-300 µL	300 µL	21 and 22	55019-09	100 µL	200 µL	21 and 22
55019-40	100 µL-1 mL	1000 µL	21 and 22	55019-11	200 µL	200 µL	21 and 22
				55019-13	250 µL	300 µL	21 and 22
				55019-15	300 µL	300 µL	21 and 22
				55019-17	500 µL	1000 µL	21 and 22
				55019-19	1 mL	1000 µL	21 and 22



SoftGrip Single Channel Pipette Carousel

Keeps pipettes safe and within easy reach. Constructed of heavy, solvent-resistant plastic.

Part No.	Description
53578-01	Holds Six Single and Multi-Channel Pipettes



Calibration Labels

Quickly record and view calibration data with labels placed on your pipettes.

Part No.	Description
1987-01	Calibration Labels, 120/pk



Calibration Key

For simple calibration in the field.

Part No.	Description
34002	Calibration Key

SoftGrip™
Pipettes

Softouch™
Electronic Pipettes

AdvanTip™
Pipette Tips

SoftAide™
Pipette Controller

SoftAide™
Quik Dispenser

Replacement Parts
and Accessories