

LINERS FOR AGILENT INSTRUMENTS

Style	Length x OD (mm)	Qty	Item No.
2 mm Splitless, Quartz	79 x 6.5	10/pk	208111-10
4 mm Split / Splitless	79 x 6.3	10/pk	208112-10
4 mm Split / Splitless Quartz	79 x 6.3	10/pk	208122-10
Split Liner, 4 mm ID with Cup	79 x 6.3	10/pk	208114-10
Single Taper, 4 mm ID with Deactivated Glass Wool	79 x 6.5	10/pk	208120-10

LINERS FOR PERKINELMER INSTRUMENTS

Style	Length	Qty	Item No.
4 mm ID Split for Autosystem	92 mm	10/pk	208232-10
4 mm ID Split for Autosystem, Deactivated Glass Wool	92 mm	10/pk	208233-10
2 mm ID Splitless for Autosystem	92 mm	10/pk	208234-10
Dimple Splitter	100 mm	10/pk	208212-10
Splitless	100 mm	10/pk	208214-10

LINERS FOR THERMO INSTRUMENTS

Style	Length	Qty	Item No.
Split, 3 mm ID	105 mm	10/pk	208611-10
Split, 5 mm ID	105 mm	10/pk	208612-10
Splitless, 3 mm ID	105 mm	10/pk	208615-10
Splitless, 5 mm ID	105 mm	10/pk	208616-10
For Wide Bore Columns	105 mm	10/pk	208617-10
Cup Liner	105 mm	10/pk	208618-10
PTV Liner, 2 mm ID	105 mm	10/pk	208620-10

LINERS FOR VARIAN/BRUKER INSTRUMENTS

Style	For Injector	Length x OD (mm)	Qty	Item No.
2 mm Splitless, Quartz	1177	79 x 6.5	10/pk	208111-10
4 mm Split/Splitless	1177	79 x 6.3	10/pk	208112-10
4 mm Split/Splitless, Quartz	1177	79 x 6.3	10/pk	208122-10
Split Liner, 4 mm ID with Cup	1177	79 x 6.3	10/pk	208114-10
Single Taper Liner, 4 mm ID with Deativated Glass Wool	1177	79 x 6.5	10/pk	208120-10

LINERS FOR SHIMADZU INSTRUMENTS

Style	Length	Qty	Item No.
Split Liner 14A	99 mm	10/pk	208414-10
Splitless Liner 14A	99 mm	10/pk	208416-10