

PEERLESS 120

PHYSICAL CHARACTERISTICS

Phase	Silica Type	Family	% C	Pore Size Å	Particle Sizes (µm)	End-Capped Yes/No	Surface Area m ² /g	USP Code
<i>Peerless 120</i>	Type B (99.999%) Ultra High Purity	C18	17%	120	3,5,10	Y	300	L-1
		C8	9%	120	3,5,10	Y	300	L-7
		AQ	15%	120	3,5	Y	300	L-1

Part No.	Range
2546571	25cm x 4.6mm x 5µm Peerless 120 C18
1546571	15cm x 4.6mm x 5µm Peerless 120 C18
1046571	10cm x 4.6mm x 5µm Peerless 120 C18

2546572	25cm x 4.6mm x 5µm Peerless 120 C8
1546572	15cm x 4.6mm x 5µm Peerless 120 C8
1046572	10cm x 4.6mm x 5µm Peerless 120 C8

2546561	25cm x 4.6mm x 5µm Peerless AQ
1546561	15cm x 4.6mm x 5µm Peerless AQ

These columns are also offered in other geometries. Please enquire with us for unlisted part numbers.