

PEERLESS BIO

PHYSICAL CHARACTERISTICS

Phase	Silica Type	Family	% C	Pore Size Å	Particle Sizes (µm)	End-Capped Yes/No	Surface Area m ² /g	USP Code
<i>Peerless Bio</i>	Type B (99.999%) Ultra High Purity	C18	8%	300	3,5,10	Y	100	L-1
		C8	5%	300	3,5,10	Y	100	L-7
		C4	3%	300	3,5,10	Y	100	L-1

Part No.	Range
2546541	25cm x 4.6mm x 5µm Peerless Bio C18
1546541	15cm x 4.6mm x 5µm Peerless Bio C18

2546542	25cm x 4.6mm x 5µm Peerless Bio C8
1546542	15cm x 4.6mm x 5µm Peerless Bio C8

2546543	25cm x 4.6mm x 5µm Bio C4
1546543	15cm x 4.6mm x 5µm Bio C4

These columns are also offered in other geometries. Please enquire with us for unlisted part numbers.